

PWSpY: A Python library dedicated to the analysis of Partial Wave Spectroscopic Microscopy data.

Nicholas M. Anthony¹ and Vadim Backman¹

¹ Department of Biomedical Engineering, Northwestern University, Evanston, IL, USA.

DOI: [10.21105/joss.03957](https://doi.org/10.21105/joss.03957)

Software

- [Review](#) ↗
- [Repository](#) ↗
- [Archive](#) ↗

Editor: [Jeff Gostick](#) ↗

Reviewers:

- [@dxm447](#)
- [@pr4deepr](#)

Submitted: 21 November 2021

Published: 03 December 2021

License

Authors of papers retain copyright and release the work under a Creative Commons Attribution 4.0 International License ([CC BY 4.0](#)).

Summary

Partial Wave Spectroscopic (PWS) Microscopy is a novel spectroscopic microscopy technique that quantitatively detects the nanoscale mass density distribution within a biological sample. PWS microscopy achieves sensitivity to nanoscale structures within biological cells by using the spectroscopic content of microscope images, and quantitatively measures nanoarchitectural changes in cells associated with carcinogenesis([Chandler et al., 2016](#)). Additionally, PWS does not require exogenous labels and thus works even in unstained samples.

PWSpY is a Python module dedicated to the analysis of PWS data. It includes a full suite of tools that prove useful in analyzing experimental data. This includes object-oriented representations of the raw data involved as well as the analyzed output data and auxiliary data used such as ROIs, calibration files, and automated imaging metadata produced by the PWS Acquisition([Anthony, 2020](#)) plugin for Micro-Manager([Edelstein et al., 2014](#)). A single PWS measurement can be treated as a 3D image cube with dimensions (x, y, and wavenumber) and a single PWS experiment may consist of thousands of such images colocalized with widefield fluorescence, confocal microscopy, or other interferometric image cubes. With PWSpY, it is trivial to skip the basics of preprocessing and get to the heart of extracting meaningful results from your experimental data. Basic operations such as normalization, camera non-linearity compensation, and calibration can all be handled with the call of a single method so you can trust that they have been performed without error. Additionally, the library provides a means for conveniently loading and storing auxiliary data such as ROIs, notes, and analysis results. Utility functionality for generating visualizations, automatic colocalization, basic modeling of reflectance, loading an manipulation of Micro-Manager position list files, and more are provided in the utility subpackage.

Statement of need

Any analysis of raw data generated by PWS microscopy or other related interferometric imaging modalities([Gladstein et al., 2019](#)) requires loading image data and metadata from a large variety of file formats and performing complex pre-processing steps before any of the real analysis begins. Up until now a great many single-use MATLAB scripts have been written to perform these tasks but a comprehensive and reusable library has never been developed. Any minor variations in how this pre-processing or analysis is performed can result in major differences in final analysis results and any script that fully implements all of the required processing becomes so long as to be unreadable. PWSpY provides an object-oriented interface for performing all common file I/O, pre-processing, and analysis tasks related to PWS. This allows users to write succinct and readable scripts/software that can be trusted to process data correctly. Automated unit testing helps ensure that as the code continues to be improved upon the analysis results will stay the same. In order to guarantee that analysis results

can always be repeated, results that are saved to file include metadata that fully describes the analysis settings, calibration files, and the Git revision SHA of PWSpy that were used in analysis.

The data processing steps described in previous publications (Lusik Cherkezyan et al., 2017), (Gladstein et al., 2019), (L. Cherkezyan et al., 2013) can be found written out clearly in Python under the `analysis` package. This library provides the backend code for the user interfaces found in `PWSpy_GUI` (Anthony, 2021). PWSpy was designed to be used by any researchers working with PWS data. It has already been used in one scientific publication (Li et al., 2021) with many more on the way. PWSpy takes care of all the low-level details of PWS analysis so that researchers can focus on the big picture.




Figure 1: Example images produced by PWSpy

Description

Analysis

A core aim of PWSpy is to make the analysis of PWS and other related interferometric measurements uniform, understandable, repeatable, and bug free. In order to make the procedures involved in each analysis pipeline in the `pwspy.analysis` package easier to understand the functionality of each type of analysis is split between a set of three classes: an analysis class that performs the actual data analysis, a settings class that defines the permissible adjustable

parameters of the analysis, and a results class which stores the output of an analysis performed on a single measurement. The outline of behavior that each of these classes must implement is defined in a set of three abstract base classes. For example, each analysis class must implement a `run` method which takes a single raw measurement as input and returns a single analysis result alongside a list of warnings indicating potential issues with the input data or analysis settings.




Figure 2: Inheritance diagram showing how the abstract base classes defining analysis pipelines, analysis settings and serializable analysis results relate to the implementations used to perform various types of image analysis

In order to help produce reproducible research each analysis result stores information about the analysis settings that were used to produce it as well as information about any calibration measurements that were used to pre-process the data and the Git revision SHA of PWSpy.

Data Types

Each type of raw experimental data that is handled by PWSpy is represented by a set of two classes, a metadata class which can provide useful information about the experimental data without loading the full image cube into memory, and a data class which facilitates operating on the 3D image data itself. In order to encourage writing code that will not be strongly coupled to the implementation details of a specific measurement type, abstract base classes for metadata and image data are defined. It is encouraged that code external to the `pws.py.dataTypes` package is written in terms of these abstract base classes when possible.

In addition to the metadata and image data classes, the `pws.py.dataTypes` package also contains classes to represent auxiliary data such as an `Roi` used to select data from a specific range within an image cube or an `Acquisition` used to group together multiple colocalized images of different modalities with their corresponding `Rois`, and any notes that have been written about the `Acquisition`.

Utility

While the core functionality of PWSpy is to facilitate the analysis of data from PWS microscopy and related interferometric imaging techniques, the `utility` package contains functionality not directly related to this task but still vital in producing well calibrated results, coregistering images from different imaging modalities, and performing higher order analysis of large-scale experiments. Below is a brief outline of two applications handled by subpackages of `pws.py.utility`, for a complete description of the functionality provided by the `utility` package please refer to the documentation.

88 **Achieving accurate spectroscopic reflectometry measurements on commercial micro-**
89 **scopes**

90 In order to accurately quantify interferometric variations in a sample's reflectance spectra
91 one must first be able to accurately quantify the reflectance of a sample. Even the highest
92 end microscope objectives will reflect some portion of the incident light back to the camera
93 when used in an epi-illumination configuration. In many cases for live-cell imaging where
94 the average sample reflectance is $\sim 0.4\%$ this unwanted internal reflectance of the microscope
95 can be even greater than the sample reflectance depending on factors like the alignment of
96 the aperture and field stops of the microscope. In order to achieve accurate reflectometry
97 measurements this internal reflectance must be accurately measured and then subtracted
98 from experimental data during pre-processing. If we model the light intensity collected at the
99 camera to be $I = I_0(R_{\text{internal}} + R_{\text{sample}})$ where I_0 is the illumination intensity and R_{sample}
100 is the sample reflectance then with measurements of two reference samples, each with a well
101 known reflectance, the internal reflectance can be calculated as $R_{\text{internal}} = \frac{R_1 I_2 - R_2 I_1}{I_1 - I_2}$

102 The `pwspy.utility.reflection.extraReflectance` subpackage provides functionality to
103 conveniently extend this simple calculation to produce an R_{internal} image cube from datasets
104 consisting of images from an arbitrary number of reference reflectances. The output of this
105 calculation is vital to the analyses contained in `pwspy.analysis`, it's data handling is managed
106 by the `pwspy.dataTypes.ERMetaData` and `pwspy.dataTypes.ExtraReflectance` classes.




Figure 3: A: Internal reflectance spectra for a single X, Y position in the microscope field of view. The final output is calculated as a weighted mean of the calculation corresponding to each pair of the available reference reflections. B: Before accounting for the internal reflectance of the microscope the ratios between measurements of two reflective materials disagree with what would be predicted by the Fresnel equations by a factor of 1.8. C: After subtracting the internal reflectance from the measurements experiment agrees well with prediction.

107 Parsing of acquisition engine metadata

108 While the core focus of PWSpy is on handling the details of accurately analyzing a single PWS
109 microscopy image cube, any real biology experiment will consist of hundreds or thousands of
110 images spanning multiple imaging conditions and cell plates. To make the acquisition of this
111 data easier the PWS Acquisition(Anthony, 2020) plugin for Micro-Manager provides flexible
112 automated imaging presented to the user as a set of possible automation steps and a drag-
113 and-drop tree which defines the order that steps will be executed in. Examples of possible
114 types of steps that could be included in an automated experiment definition are: a time series
115 of sub-steps, execution of substeps at a pre-defined set of XYZ coordinates, acquisition of a
116 PWS image, a pause to wait for user input, etc.

117 In order to minimize the amount of custom code that a data analyst must write for each new
118 experiment in order to locate the data of an individual analyzed image in the larger context of
119 the entire experiment, the `pws.py.utility.acquisition` subpackage contains functionality
120 for parsing the metadata and log files saved by the acquisition engine. An instance of the
121 `SequenceAcquisition` class binds each individual `pws.py.dataTypes.Acquisition` to
122 a `SequencerCoordinate` which defines a unique location in the experiment sequence of
123 automation steps, each of which is represented by an instance of `SequencerStep`.

124 Availability

125 PWSpy is free and open source. It is published under the GNU General Public License v3. You
126 can download the source code at <https://github.com/BackmanLab/PWSpy>. Online docu-
127 mentation and examples can be found at <https://pws.py.readthedocs.io>. PWSpy can be in-
128 stalled using Pip or using Conda from the “backmanlab” Anaconda Cloud channel.

129 Acknowledgements

130 The authors would like to thank many of the early users of PWSpy for providing feedback and
131 bug reports as the project has developed.

132 References

- 133 Anthony, N. M. (2020). PWS acquisition: A micro-manager plugin for the automated
134 acquisition of multimodal interference-based imaging. In *GitHub repository*. GitHub.
135 <https://github.com/nanthony21/PWSAcquisition>
- 136 Anthony, N. M. (2021). PWSpy GUI: A collection of graphical user interfaces for the analysis
137 of multimodal interference-based image data. In *GitHub repository*. GitHub. https://github.com/nanthony21/pws.py_gui
- 138
- 139 Chandler, J. E., Cherkezyan, L., Subramanian, H., & Backman, V. (2016). Nanoscale refrac-
140 tive index fluctuations detected via sparse spectral microscopy. *Biomedical Optics Express*,
141 7(3), 883. <https://doi.org/10.1364/boe.7.000883>
- 142 Cherkezyan, L., Capoglu, I., Subramanian, H., Rogers, J. D., Damania, D., Taflove, A.,
143 & Backman, V. (2013). Interferometric spectroscopy of scattered light can quantify the
144 statistics of subdiffractional refractive-index fluctuations. *Physical Review Letters*, 111(3).
145 <https://doi.org/10.1103/physrevlett.111.033903>
- 146 Cherkezyan, Lusik, Zhang, D., Subramanian, H., Capoglu, I., Taflove, A., & Backman, V.
147 (2017). Review of interferometric spectroscopy of scattered light for the quantification of
148 subdiffractional structure of biomaterials. *Journal of Biomedical Optics*, 22(3), 030901.
149 <https://doi.org/10.1117/1.jbo.22.3.030901>
- 150 Edelstein, A. D., Tsuchida, M. A., Amodaj, N., Pinkard, H., Vale, R. D., & Stuurman, N.
151 (2014). Advanced methods of microscope control using uManager software. *Journal of*
152 *Biological Methods*, 1(2), e10. <https://doi.org/10.14440/jbm.2014.36>
- 153 Gladstein, S., Almassalha, L. M., Cherkezyan, L., Chandler, J. E., Eshein, A., Eid, A.,
154 Zhang, D., Wu, W., Bauer, G. M., Stephens, A. D., & al., et. (2019). Multimodal
155 interference-based imaging of nanoscale structure and macromolecular motion uncovers
156 UV induced cellular paroxysm. *Nature Communications*, 10(1). <https://doi.org/10.1038/s41467-019-09717-6>
- 157

158 Li, Y., Eshein, A., Virk, R. K. A., Eid, A., Wu, W., Frederick, J., VanDerway, D., Gladstein,
159 S., Huang, K., & Shim, A. R. et al. (2021). Nanoscale chromatin imaging and analysis
160 platform bridges 4D chromatin organization with molecular function. *Science Advances*,
161 7(1). <https://doi.org/10.1126/sciadv.abe4310>

DRAFT